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Melting Point Meter M5000

User Manual

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Index

English

1. Description M5000	19
2. Safety Precautions	20
3. Unpacking and Assembly	20
4. Melting point analysis	21
5. Structure process / operator programme	22
6. Starting the melting point analysis	23
7. Reaching the melting point	25
8. Subsequent retrieval of the last melting point	26
9. Setting of date, time and language	27
10. Technical data	28
11. Maintenance	29
12. Recovery and Recycling	29
Warranty conditions	30

1. Description M5000

2. Safety precautions

- Do not cover ventilation slots.
- Do not insert any objects into ventilation slots.
- No humidity may intrude into the device.
- During analysing the melting point, the measuring head heats up extremely around the ventilation slots. Do not touch the measuring head in this area.
- Only use suitable capillary tubes with an external diameter of 1.4 mm for the determination of the melting point.
- When working with chemicals, follow safety precautions.
- When maintaining the device, e.g. changing the fuse, cleaning with a dry or any slightly wet towel, always make sure that the device is disconnected from electricity.
- Only authorized personnel or the manufacturer may carry out servicing and maintenance.

3. Unpacking and assembly

- Carefully unpack the device.
- Immediately report any transport damage occurred.

Find a suitable place for the device, i.e. avoid direct sunlight and vicinity of heating devices. Connect power plug to 230 V / 50 Hz alternating current. Before initial measuring the device should be switched on for about 15 min to reach ideal working temperature.

4. Melting point analysis

With the device you will be able to determine the melting point of powdery substances up to 400 °C.

- ▶ Switch on the mains switch on the back of the device.
You will see the current oven temperature on the display and you will hear a signal tone.
- ▶ By pressing the arrow key, select the desired starting temperature (approx. 3 °C below expected melting point).
- ▶ When having reached the selected starting temperature, carefully insert a filled capillary tube (filled with approx. 4 mm) with the opening facing upwards into the drill hole of the guiding head until you feel the capillary tube touching the end. After having reached the melting point, the device falls back to the selected temperature and you will also hear a signal tone. The melting point is displayed and can be printed via the printing key.
- ▶ During cooling down to the starting temperature “auto/waiting” will flash up on the display. After the “waiting” display has disappeared, the same measuring can be repeated.
- ▶ If there is no new measuring within 5 min, the device automatically cools down (fan) to 30 °C.
- ▶ Do not switch off the device before reaching 30 °C, because then the fan cooling will also be switched off.
- ▶ If you wish to carry out a measuring while the cooling is operating, you can directly select the temperature using the arrow keys. The cooling will switch off automatically.
- ▶ If the desired starting temperature is below the current temperature, you have to wait until the cooling has reduced the temperature.
- ▶ Do not attach protective cover before the device has cooled down.

5. Structure process / operator programme

5.1 Starting display

actual temperature
25.8 °C

The display always appears when switching on the device.

If the temperature is above 35 °C, the fan will start automatically.

If the temperature is below 30 °C, the fan will switch off.

If one of the arrow-keys is pressed, you continue with 5.2 *Selection of the target temperature.*

5.2 Selection of the target temperature

select temperature
120.0 °C

Adjustment takes place in steps of 1 °C.

The menu can only be left by turning the device off and on or by time-out.

After 5 min (without pressing any key) time-out appears and the device will be restarted (see 5.1 Starting display).

When choosing a temperature, the fan will always be switched off.

If a starting temperature higher than the current oven temperature is selected, inserting the capillary tube will start the normal process of melting point analysis.

If the starting temperature is lower than the current oven temperature, an error message will be displayed for 3 s when the capillary tube is inserted.

Selected temperature
> oven temperature

After that, the device will return to the starting display and the procedure will be cancelled. The fan will be started according to the logic.

The latest selection of temperature can only be erased when turning off the device.

6. Starting the melting point analysis

6.1 In principle

The starting of a melting point analysis will only be effected by inserting a filled capillary tube in the oven, if a starting temperature has previously been selected and reached.

In the case of an unfilled capillary tube, an error warning will be displayed. It will be automatically deleted after the empty capillary tube has been removed.

empty capillary tube

After inserting, measurement will automatically start. When the capillary tube is inserted, the starting menu is automatically shown.

If there is already a capillary tube in the oven before the starting temperature has been selected, a message on the display will ask you to remove the capillary tube.

Please remove
capillary tube

Only then (automatically after removing the capillary tube) the selecting menu will be shown.

The last selected starting-temperature remains in memory until the device is turned off.

The last melting point determined will be saved until the device is turned off or a new melting point analysis is started.

6.2 Express heating-up

express heating
Temp. 78.8 °C

After selecting the starting temperature, express heating up will automatically be started after pressing the start key. The current oven temperature will be displayed.

After reaching the selected starting temperature (the device decreases heating-up speed approx. 3 °C before) a signal tone will be heard for about 5 seconds.

Express heating up can be cancelled by pressing one of the arrow keys. The selection of the target temperature (see chapter 5.2) is started with the arrow keys. The fan will not be activated! The printing key has no function.

6.3 Slow heating-up

1 °C/min → 2 %
Temp. 121.3 °C

Slow heating-up means approx. 1 °C per minute. The process is automatically started by inserting a filled capillary tube when the selected starting temperature has been reached. The melting process is displayed in °C and % from the beginning of melting until the melting-point is reached. The current target temperature of the melting-oven is displayed.

In case the melting point is not reached within 30 min, the slow heating-up process will be interrupted by time-out and an internal error value will be saved as melting point.

By pressing the arrow keys, the process can be interrupted at any time. The device will then jump to menu item 5.2 *selection of target temperature*. No fan will be activated! The printing key has no function.

7. Reaching the melting point

melting point
121.9 °C auto

A signal tone will sound three times when the melting point has been reached.

The determined melting point will be shown on the display. The device will automatically set the melting oven back to the selected starting temperature. During this time the term "auto" will flash up in the lower right corner of the display.

When the starting temperature has been reached once more, the signal tone will sound and the term "auto" will remain in the lower right corner of the display.

If a new capillary tube is inserted during the time the term "auto" is flashing up, no error message will be displayed. As soon as the starting temperature has been reached, the next melting point analysis process will immediately be started.

If a new melting sample is inserted within 5 min, the measuring process will be repeated.

If one of the two arrow keys is pressed within 5 min, the selection of the target temperature will be activated (see chapter 5.2) No fan will be turned on!

If the printing key is pressed within 5 min, you can print on the optionally attached printer. This process can be repeated as often as you want.

If no keys are pressed for about 5 min, a time-out will occur and the device will automatically jump to the starting display. The last melting point remains in memory and can be displayed at a later time by pressing the printing key.

8. Subsequent retrieval of the last melting point

actual temperature
25.8 °C

If you missed the time for reading the melting point and the device has returned to the starting display, the last melting point can be displayed once more by pressing the printing key.

melting point
121.9 °C

The display will appear for about 5 sec. After that the device will return to the starting display. If the printing key is activated once more within the 5 s, the print-out will be created with the attached printer (if available). This can be repeated as often as you want.

The term "auto" will not be shown on the display, since no melting point analysis can be started in this case inserting a capillary tube does not have any effect.

If an internal error value is determined as melting point, the melting point analysis was prematurely aborted by pressing the arrow keys (clear text display of the error for 5 s).

The melting point could not be determined due to a time-out (10 min.), i.e. the melting point could be 10 °C higher than the selected starting temperature (clear text display of the error for 5 s).

abort,
no melting point

9. Setting of date, time and language

After pressing the time-set-key on the right hand side of the display, the set-up menu appears. The settings will be made in the following order:

min
h
weekday
year
month
day
language
print

return with:

english12:38
Mo 31.01.2009

A flashing cursor will appear on the display. With each activation of the time- / set-key the cursor will jump to the next item. The value can be changed resp. set with the arrow keys. When setting the date you should start with year, then month and at least day and week-day. When the month is set, the day always shows 1. Hour and minute are set separately.

When all settings are completed, the menu is left and all settings are saved by pressing the printing key. The settings will also remain after turning off the device.

The device includes a battery. This is why the date and time are also updated in the off mode for about 1 year without turning on the device.

Calibration parameters as well as language are stored in EEPROM and therefore remain for a longer period.

10. Technical data

Power supply	230 V
Power consumption	40 VA
Safety classified	IP 20
Measuring range	25–400 °C
Accuracy	25–200 °C, ± 0.3 °C
	200 °C–400 °C, ± 0.5 %

Heating-up rate during measuring 1 °C/min.
It is recommended to have the device calibrated every two years.

11. Maintenance

Like all precision devices, your Melting point meter should be treated with care. Protect it against dust, condensation, vibration and direct solar radiation. Please do not apply acid detergents to the instrument!

Important information!
Corrosive chemicals can harm the varnish or plastic housing!

Important information!
In case of defect, the device must be repaired by the company A.KRÜSS or an authorised workshop.

12. Recovery and Recycling

The device can be an important source of raw materials. Please do not dispose of as waste, but collect separately for the recycling and recovery of the contained materials. If disposed improperly, the materials may be damaging to the environment and human health.

The manufacturer of the device, A.KRÜSS Optronik GmbH, collects, uses and recycles the contained raw materials. However this recovery requires your support.

If you decide to dispose of this melting point meter, please do not try to open it up or to use parts of it in any other way than described in this manual, but return the device to the dealer you purchased it from.

The dealer should take the device back free of cost.

The recovery of the raw materials is effected with respect to the European guideline 2002/96/EC and any other applicable guideline.

Warranty conditions

A.KRÜSS Optronic guarantees the materials and workmanship of KRÜSS instruments for a period of 24 months from the date of shipping. A.KRÜSS Optronic will repair or replace defective devices within this period, given that they fall under terms of the guarantee. The device must be sent back to A.KRÜSS Optronic for warranty repairs or service. Shipping from the customer is at the expense of A.KRÜSS Optronic for warranty repairs. In all other cases, expenses must be paid by the customer.

A.KRÜSS Optronic guarantees that the hardware specified by A.KRÜSS Optronic for use with this device will function without error if used according to the instructions provided in our manual.

A.KRÜSS Optronic does not guarantee error-free and uninterrupted operation of the device or the accuracy of this instructions manual. A.KRÜSS Optronic is not liable for consequential damage.

Warranty limitations:

This warranty does not cover errors and damage arising due to improper handling, through modification, misuse, operation above and beyond the specified environmental conditions or through unauthorised maintenance.

Further claims will not be accepted or recognised. A.KRÜSS Optronic expressly provides no guarantee for the workability or economic use in specific application cases.

A.KRÜSS Optronic reserves the right to make changes to this instruction manual and the technical data of the device in question at any point in time.

All devices may only be transported when properly packed in the complete original packaging. Request replacement packaging from your supplier if necessary.

A.KRÜSS Optronic GmbH
Alsterdorfer Straße 276–278
22297 Hamburg | Germany
Tel +49-(0)40-51 43 17-0
Fax +49-(0)40-51 43 17-60
E-Mail info@kruess.com
Web www.kruess.com

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Tel +49-(0)40-51 43 17-0
Fax +49-(0)40-51 43 17-60

E-Mail info@kruess.com
Web www.kruess.com